

Login / Register

Contact

I'm looking for

0 Items \$0.00

Home Prosthetics Orthotics Mobility Services Education




Home Orthotics Bracing & Supports Spine LO: L1-S1 **Low Profile LSO**

## Low Profile LSO

Item #: 50R151N

Classic LO (L1-S1)

Log in to order



Feedback

Product Details

Technical Specifications

**LO offers compression and stabilization from L1 to S1 with a 5:1 mechanical advantage.**

### Features

- Unsurpassed compression and support with a 5:1 mechanical advantage pulley system
- One-hand pull system compresses and cones to individual patient anatomy
- Low-profile corset and panels for increased patient compliance
- Modular step-down system is ideal for treatment pre-op, post-op, and through rehabilitation
- Mesh corset for increased patient comfort

### Indications

- Indicated for low back pain sprains and strains, lumbar disc displacement, osteoporosis, disc degeneration, disc herniation, post-operative fusion, laminectomy, discectomy, IDET procedures, spondylolisthesis and spinal stenosis.

Prosthetics

Orthotics

Mobility

Connect

