

MIAMI J® CERVICAL COLLAR

Known as the superior c-spine immobilizer on the market through multiple independent studies, the Miami J Collar's patented design meets the specialized needs of all those afflicted with cervical injuries. The Miami J offers features that enhance compliance and help speed recovery while ensuring immobilization, inhibition of skin breakdown and patient comfort.

- Sorbatex[™] padding is antibacterial and clincally proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- Bioengineered to minimize pressure points in key known areas; chin, occiput, trapezious and clavicle
- · Sizing based on 6 phenotypes to assure proper fit and optimum immobilization
- Comfort adjustment buttons for occiput and mandible, removable sternal pad for supine and swallowing

Collar Sets & Accessories (see Phenotype-driven sizing system below)

C CL L/I/ L L: 4
Super Short/Kyphotic*
Stout
XS
Tall
Regular
Short
Bariatric

Front & Back Panels (see Phenotype-driven sizing system below) **MIAMI | FRONT MIAMII BACK** SI7F

, , , , , , , , , , , , , , , , , , , ,	, , , , , , , , , , , , , , , , , , , ,	
MJ-200SFRT	MJ-400BK	Super Short/Kyphotic*
MJ-200LFRT	MJ-200LBK	Stout
MJ-250FRT	MJ-250BK	XS
MJ-500FRT	MJ-500BK	Tall
MJ-400FRT	MJ-400BK	Regular
MJ-300FRT	MJ-300BK	Short

Miami J Front and Back Panels may be interchanged to accommodate patient anatomy, e.g., a Small Front may be used with a Large Back.

Choose collar size based on patient neck height, not patient body size.

Standard sizes: 60% of the population

PHENOTYPE SIZING CHART

Specialty sizes: 30% of the population



ankylosing spondylitic









- Young women
 Adolescents



- mature women

 Thin, mature male
- Short 30% (Shrt) All other adults? Standard size for Short-necked women

Phenotype-driven sizing system makes size selection simple while drastically reducing margin of error in fitting a patient. Each collar has a bioengineered front with a specific back designed to work together for superior immobilization.

^{*} Super Short/Kyphotic collars are packaged with Medium Back Panels.