



NEW OPTIONS SPORTS

**ORTHOPEDIC AND SPORTS
MEDICINE CATALOG**



Distributed by
Frontline Medical Services



www.newoptionssports.com

Phone 800.872.5488 • Fax 800.455.5488

INTRODUCTION AND INDEX



New Options Sports has been in the orthopedic and prosthetic arena for over 40 years. Located in Farmers Branch, Texas, we proudly manufacture American-made orthopedic sports medicine bracing as well as prosthetic limb covers and suspension sleeves for amputees.

In addition to products for adults, we offer pediatric products. John Scott, President and owner of New Options Sports has long supported the idea that, "Kids need braces too!"

For both adult and pediatric products, New Options Sports is one of the few orthopedic manufacturers offering a custom service program for hard to fit patients. Turn around time for a custom product is 3- 4 business days and additional material colors are available.

Quality and design of New Options Sports products are second to none. Only the best latex free materials are utilized. At New Options Sports, we understand the benefits of neoprene. Neoprene uses the body's natural warmth to provide concentrated heat to the affected area. Our 4-way stretch neoprene provides superior compression and support.

We also manufacture a line of identically designed products made of Koolflex™ material. Koolflex™ is a light, breathable knit fabric which wicks perspiration away from the body. Koolflex™ is an alternative product for those who have to wear supports for long periods of time or are in warm climates.

New Options Sports distributes through specialty distributors, orthotic and prosthetic distributors, retail outlets, DME facilities and national athletic supply companies. We are a proud supplier to military facilities and veterans hospitals throughout the United States. We also manufacture products for internet sales and global distribution companies. Our custom logo program is always available for customers wishing to promote their brand to their clients and patients.

MADE in the USA is a very important aspect of our company. We manufacture our products in the United States. American pride, integrity and a wealth of knowledge have made New Options Sports what it is today.

INDEX

OSTEOARTHRITIS BRACING	1	ANKLE SUPPORTS	8
HINGED KNEE SUPPORTS	2-3	SPORTS MEDICAL PRODUCTS	8-9
KNEE SLEEVES	4	LUMBAR SUPPORTS & KOOLFLEX™, STAYS, HINGES	10
PATELLA STABILIZERS	4-5	SUSPENSION SLEEVES and C-ME™ COVERS	11
WRIST SUPPORTS	6	OA PART NUMBERS PULL-UP AND WRAPAROUND	12
ELBOW SUPPORTS	7	FITTING GUIDE	13



Proudly made in the USA

OSTEOARTHRITIS BRACING

OA and OA WRAP

RIGHT HINGE FITS:
RIGHT MEDIAL
LEFT LATERAL

LEFT HINGE FITS:
LEFT MEDIAL
RIGHT LATERAL




L1843 Recommended

“NEW OPTION” OA BRACE

Features:

- Compressive buttressed knee sleeve, assuring correct joint alignment and proper donning
- Soft construction for stability without “over bracing”
- Simple adjustment and fit with proximal and distal closures
- Strut paddles for rotation control and load distribution
- Pre-set 7° offset for optimal counter force and unloading
- Optional 10° & 13° condyle pads for additional counterforce, if needed
- Push technology unloads joint line opposite of hinge side
- Low profile, lightweight design for discreet support beneath clothing
- Internal sharkskin to control migration



L1843 Recommended

OA WRAP AROUND BRACE

Features:

- Wrap-around, front closure style for easy donning.
- Sharkskin patella buttress assists with patella tracking
- Same hinge features as the original NOS OA brace above
- Brace made with **KoolFlex** material
- Internal sharkskin to control migration



Removeable Patella Buttress

See page 12 for sizing chart and a complete list of part numbers



Proudly made in the USA

HINGED KNEE SUPPORTS

KC64-NOS KNEE MATE™



L1832 Approved

KNEE MATE™ Wraparound

Indications:

- Medial and lateral collateral instability
- Anterior and posterior cruciate instability
- Post-operative management
- Arthritis

Features:

- Breathable KoolFlex material
- Patella opening for positioning and comfort
- Removeable upper & lower circumferential straps with elastic strap segment for improved fit and comfort
- Large, comfortable popliteal opening
- Stainless steel hinges w/flexion and extension stop kit
- Special cordura stress covers over center, proximal and distal hinge pocket for reinforcement

Product Variation:

KC64-SL **KoolFlex** SuperLite with simple hinge -L1810
 K64-PC in 3/16" neoprene with covered PolyCentric hinges
 KC64-PC in **KoolFlex**
 K64-MP in 3/16" neoprene w/ covered Multi-Position hinges
 KC64-MP w/ covered Multi-Position hinges in **KoolFlex**
 K64-HT in 3/16" neoprene w/ high tech, ROM hinges
 KC64-HT w/ high tech, ROM hinges in **KoolFlex**
 K64-NOS in 3/16" neoprene w/ NOS hinges
 KC64-NOS w/ NOS hinges in **KoolFlex**
 K64-E in 1/8" neoprene with spiral stays
 K65-HT 17" brace in 3/16" neoprene w/ high tech, ROM hinges
 KC65-HT 17" brace w/ high tech, ROM hinges in **KoolFlex**
 K65-NOS 17" brace in 3/16" neoprene w/ NOS hinges
 KC65-NOS 17" brace w/ NOS hinges in **KoolFlex**

KC65-NOS KNEE MATE™



L1832 Approved

KC67-NOS THE HYBRID™



L1832 Approved

"THE HYBRID™"

Features:

- **KoolFlex™** material open-cell technology allows patient to stay cool in brace
- Patella opening for positioning and comfort
- Anterior wrap around closure on top, closed bottom around calf area - better fit for the cone-shaped leg
- Removable upper & lower circumferential straps with elastic strap segment for improved fit and comfort
- Open popliteal
- Adjustable and covered stainless steel hinges with flexion and extension stop kit
- Fits right or left knee

Product Variations:

K67-PC in 3/16" neoprene with PolyCentric hinge
 KC67-PC in **KoolFlex**
 K67-MP in 3/16" neoprene w/ covered multi-position hinges
 KC67-MP w/ covered multi-position hinges in **KoolFlex**
 K67-HT in 3/16" neoprene w/ high tech, ROM hinges
 KC67-HT w/ high tech, ROM hinges in **KoolFlex**

KC68-HT- 17" Hybrid w/high tech, ROM hinges in **KoolFlex**



U.S. Patent 7217249

KC68-NOS THE HYBRID™



L1832 Approved



HINGED KNEE SUPPORTS



K1-U-MP PULL-ON COVERED HINGED KNEE

Features:

- 3/16" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Universal felt buttress
- Upper & lower circumferential straps
- Additional straps above and below patella for support
- Covered Multi-Position hinges
- Boot pulls
- Closed popliteal
- **L1832 approved**



K2-U-MP PULL-ON PADDED HINGED KNEE

Features:

- 3/16" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Universal felt buttress
- Upper & lower circumferential straps
- Covered and padded Multi-Position hinges
- Boot pulls
- Closed popliteal
- **L1832 approved**

Product Variations:

KC2-U-MP in KoolFlex



K12-NOS Over and Under™

Features:

- Breathable KoolFlex material
- Anterior closures for easy donning and doffing
- Fits universal left or right
- Comfortable large popliteal opening
- Adjustable universal patella buttress
- Removable upper and lower circumferential straps w elastic strap segment
- Range of motion hinges w/flexion and extension stops
- **L1832 approved**

Product Variations:

K12-NOS - with NOS Hinges in 4-way stretch neoprene



K34-PC PULL-ON "SUPERLITE™" KNEE BRACE

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Upper & lower circumferential straps
- Polycentric aluminum, low-profile hinges
- Boot pulls
- Open popliteal
- Large size weighs 8.6 oz
- **L1820 approved**



K42-HT WRAPAROUND REHAB KNEE BRACE WITH ROM HINGES

Features:

- 1/8" nylon, 4-way stretch neoprene
- Patella opening for positioning and comfort
- Popliteal opening
- Two upper and lower circumferential straps
- Flip-top ROM hinge for locking out or set for flexion and extension control
- Boot pulls



- **L1832 approved**



K45-MP PULL-ON SWEDISH KNEE BRACE

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Universal neoprene crescent buttress
- Upper & lower circumferential straps
- Additional straps above and below patella for support
- Covered Multi-Position hinges
- Open popliteal
- **L1832 approved**



Proudly made in the USA

KNEE SLEEVES & SUPPORTS

K9-OP KNEE SUPPORT

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning & comfort
- Anterior oval pad

Product Variation:

K9-CP w/ closed patella & oval anterior pad

K9-O w/ patella opening, no pad

KC9-O w/ patella opening, no pad in

KoolFlex

K9-C w/ closed patella, no pad

K9-OD w/open patella and 1/4" enclosed tubular donut

K9-OST open patella w/ superior tubular buttress



K10 PERFORATED UNDERSLEEVE FOR FUNCTIONAL KNEE BRACE

Features:


- 1/16" nylon perforated neoprene
- 19" Length
- Patella opening for positioning and comfort
- Provides suspension for functional knee brace
- Protects user from contact with outer brace
- Provides warmth and equal compression
- **L2397 approved**

K18 OSGOOD SCHLATTER SUPPORT

Features:

- 1/8" nylon 4-way stretch neoprene
- 1/4" v-shaped, shock absorbing pad to protect tibial tuberosity
- Upper & lower circumferential straps
- Lower strap over tibial tuberosity to apply desired pressure

Patent 362,305



PATELLA STABILIZERS

K4-F PULL-ON PATELLA FAN KNEE STABILIZER DIAMOND X THEORY

Features:

- 3/16" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Universal felt buttress
- Upper & lower circumferential straps
- Medial and lateral spiral stays
- Four 2" wide criss-cross elastic straps for added stability
- Closed popliteal



K4-U 3/16" PULL-ON PATELLA KNEE STABILIZER WITH UNIVERSAL FELT BUTTRESS

Features:

- 3/16" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Universal felt buttress
- Upper & lower circumferential straps
- Medial and lateral spiral stays
- Closed popliteal



K8-U 1/8" PULL-ON PATELLA STABILIZER WITH UNIVERSAL FELT BUTTRESS

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Universal felt buttress
- Upper & lower circumferential straps
- Medial and lateral spiral stays
- Closed popliteal



K14-D PULL-ON PATELLA STABILIZER W/ ADJUSTABLE TUBULAR BUTTRESS DONUT

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- 1/4" enclosed tubular donut
- Upper & lower circumferential straps
- Medial and lateral spiral stays
- Closed popliteal



PATELLA STABILIZERS



K15 DYNAMIC PULL-ON PATELLA STABILIZER WITH J BUTTRESS

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Patella pull-panel with adjustable straps
- 1/4" tubular, enclosed j-pad buttress
- Medial and lateral spiral stays
- Closed popliteal
- Specify right or left



K17-PC HINGED PULL-ON DYNAMIC PATELLA STABILIZER WITH J BUTTRESS

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Patella pull-panel with adjustable straps
- 1/4" tubular, enclosed j-pad buttress
- Polycentric Hinges
- Upper and lower circumferential straps
- Closed popliteal
- Specify right or left
- **L1820 Approved**



KC44-NOS DYNAMIC PATELLA STABILIZER W/ UNIVERSAL SHARKSKIN PULL BUTTRESS

Features:

- Breathable KoolFlex material
- Patella opening for positioning and comfort
- Large, popliteal opening
- Universal sharkskin buttress with dynamic pull-tabs to adjust pressure from outside of support
- Removable circumferential straps w/elastic strap segment
- Range of motion hinges with flexion and extension stops
- **L1832 Approved**

Product Variation:

K44-NOS in 1/8" neoprene w/NOS hinges
Patent 7,862,528




K46 PULL-ON PATELLA KNEE SUPPORT W/ EXCLUSIVE POSITIVE CONTROL DISTAL STRAP

Features:

- 1/8" nylon 4-way stretch neoprene
- Patella opening for positioning and comfort
- Inverted horseshoe has sewn in loop for dynamic pull
- Lower strap has guide for positive positioning and control of patella strapping
- Upper & lower circumferential straps
- Medial and lateral spiral stays
- Closed popliteal

Product Variation:

KC46 in KoolFlex



K197 SHARKTRAK™ HINGED WRAP AROUND PATELLA STABILIZER

Features:

- 1/8" UBL 1-side, nylon 1-side neoprene
- Dual Dynamic Pulls
- Light duty hinges
- Patella opening for positioning and comfort
- Two Closures w/ anchor tabs for easy application
- Universal Left or Right dynamic pull tab buttress for desired pressure on the patella
- Shark Skin buttress material grips the skin to control patella
- **L1820 Approved**

Product Variation:

K196 with spiral stays

Patent 7,959,59



K199 THE SPOONER

Features:

- Adjustable alligator clips
- 1/8" nylon 4-way stretch neoprene
- Universal Left or Right dynamic pull tab buttress for desired pressure on the patella
- Shark Skin buttress material grips the skin to control patella
- Prevents migration of functional brace
- Can be worn alone, directly on knee



Proudly made in the USA

WRIST SUPPORTS

W1 ACTION WRIST SUPPORT

Features:

- 1/8" nylon 4-way stretch neoprene
- Hook and loop closures
- Circumferential strap to limit motion
- Removable palmar splint for washing
- 6" in length
- **L3908 Approved**
- *Specify right or left hand*



W11 THUMB SPLINT

Features:

- 1/8" nylon 4-way stretch neoprene
- Hook and loop closures
- Adjustable thumb strap
- Universal size
- Fits right or left hand
- 6" in length
- **L3807 Approved**



W20 UNIVERSAL UBL WRISTLET

Features:

- 1/8" 4-way stretch neoprene, outside is UBL, inside is nylon
- Hook and loop closures
- Adjustable thumb strap
- Circumferential strap
- Removable palmar splint for washing
- 6" in length
- **L3908 Approved**
- *Specify right or left hand*



Product Variation:

W28 - 8" in length

W47 WRIST / THUMB SUPPORT

Features:

- 1/8" nylon 4-way stretch neoprene
- Hook and loop closures
- Abducted, adjustable thumb strap
- Circumferential strap
- Removable palmar splint for washing
- **L3807 Approved**
- 6" in length
- *Specify right or left hand*



W57 THUMB AND WRIST SUPPORT

Features:

- 1/8" nylon 4-way stretch neoprene
- Hook and loop closures
- *Specify right or left hand*



W56 PULL-ON THUMB AND WRIST SUPPORT

Features:

- 1/8" nylon 4-way stretch neoprene
- *Specify right or left hand*



W59 VELPLUSH™ UNIVERSAL BOOMERANG WRIST WRAP

Features:

- 1/8" VeLPlush™ neoprene
- Hook and loop closure
- Fits right or left hand
- Universal size



KOOLFLEX COMFORT WRIST SUPPORTS

WC32 - 8" KOOLFLEX™ WRIST SUPPORT

Features of both:

- 1/8" **Koolflex** breathable, hypoallergenic material
- Circumferential, perforated elastic compression strap secures palmar & dorsal splints to limit motion and ensure proper fit
- Removable dorsal and palmar splints for washing
- **L3908 Approved**
- *Specify right or left hand*

Indications:

Sprains
Strains
Carpal tunnel syndrome
Post casting
Rehabilitation
Preventative



WC33 - 10" KOOLFLEX™ WRIST & FOREARM SUPPORT



ELBOW SUPPORTS



E1 ELBOW SLEEVE

Features:

- 1/8" nylon 4-way stretch neoprene
- 10" in length
- Provides warmth, compression and protection

Product Variation:

E8 w/ covered shock absorbing pad - see Sports Medical Products



E4 TENNIS ELBOW STRAP W/ FELT PAD

Features:

- 1/8" nylon 4-way stretch neoprene
- Covered 1/4" felt pad
- Circumferential loop and lock strap w/ elastic segment



E5 TENNIS ELBOW STRAP W/ VELFOAM PAD

Features:

- 1/8" nylon 4-way stretch neoprene
- Covered velfoam pad
- Circumferential loop and lock strap w/ elastic segment



E6-U UNIVERSAL TENNIS ELBOW STRAP W/ GEL PAD

Features:

- 1/8" nylon 4-way stretch neoprene
- Covered gel pad
- Circumferential loop and lock strap w/ elastic segment
- Universal size



E9-MP HINGED ELBOW

Features:

- 1/8" nylon 4-way stretch neoprene
- Covered Multi-Position hinges
- Upper and lower circumferential straps
- L3710 Approved



E12-PC HYPEREXTENSION HINGED ELBOW

Features:

- 1/8" nylon 4-way stretch neoprene
- Covered, stainless steel polycentric hinges with hyperextension stop
- Two loop and lock criss-cross straps to help prevent hyperextension
- L3710 Approved



E15-PC HYPEREXTENSION POLYCENTRIC HINGED ELBOW W/ "E-Z-E" BICEP CLOSURE

Features:

- 1/8" nylon 4-way stretch neoprene
- Covered Polycentric hinge with hyperextension stops
- "E-Z-E" bicep closure
- One continuous bottom strap with criss-cross lock for desired angle stop
- L3710 Approved

Product Variation:

EC15-PC in Koolflex



EC16 CUBITAL TUNNEL ELBOW SUPPORT

Features:

- Breathable Koolflex material
- Power knit mesh inner sleeve
- Hook and loop closure
- Removable stays for washing
- Can be worn overnight
- L3762 Approved

Product Variation:

E16 in neoprene

U.S. Patent 7217249



Proudly made in the USA

ANKLE SUPPORTS



A30 WOOTEN-UHL (5-IN-1) FOOT AND ANKLE ORTHOSIS

Indications:

Designed to help relieve heel pain, the orthotic device applies counterforce pressure on the Achilles Tendon

- **Achilles Tendon** - place pad in superior position
- **Haglund's Disease** - for pump bump deformity, place pad in inferior position
- **Sever's Disease** - traction epiphysitis of the calcaneus, place pad in inferior position
- **Heel Spur** - invert neoprene product and place pad beneath heel w/ horseshoe pad open towards toe
- **Plantar Fasciitis** - invert neoprene product and place pad appropriately between heel and arch

Features:

- 1/8" nylon 4-way stretch neoprene
- Hook and loop ankle closure
- Removable 1/4" PPT horseshoe

U.S. Patent # 4841957



A50 COOPER II ANKLE SUPPORT

Indications:

- Lateral and medial instability
- Long-term management
- Post-operative management

Features:

- 1/16" perforated nylon 4-way stretch neoprene
- Fits right or left ankle and easily into any shoe
- Hook and loop closures
- Medial and lateral non-stretch stabilization straps
- Two stabilizing top straps
- U-shaped 15" spiral stay prevents migration
- Protective support
- Helps resist inversion and eversion
- Allows for full dorsal and plantar flexion

The Cooper II Ankle Stabilizer is bilateral and fits either left or right foot. This reduces inventory by eliminating the need for stocking lefts and rights. The medial and lateral cinch straps originate from the forefoot to better control inversion and eversion and other instabilities of the ankle. The Cooper II ankle stabilizer brace is comfortable to wear and simple to apply. The Cooper II Ankle Stabilizer brace takes the patient from the emergency room back to the field of competition.

U. S. Patent # 505062

SPORTS MEDICAL PRODUCTS



C1 CALF SUPPORT

Features:

- 1/8" nylon 4-way stretch neoprene
- 11" sleeve



E8 PADDED ELBOW SLEEVE

Features:

- 1/8" nylon 4-way stretch neoprene
- Covered, shock-absorbing pad



DC2 ICING DIGIT COVERS

Features:

- 1/8" nylon 4-way stretch neoprene
- 1 pair
- Universal size



F1 CRYO BOOTIE

Features:

- 1/8" nylon 4-way stretch neoprene
- Non-slip neoprene sole
- Finger loop for easy application
- 1 each





G1 GROIN STRAP UNIVERSAL

Features:

- 3/16" UBL 1-side, non-slip 1-side neoprene
- 6" wide contoured band
- Universal size
- Hook closures for adjustable compression



S5 THE ARM-ADILLO™ SHOULDER STABILIZER

Indications:

- Designed to limit abduction of and support shoulder girdle
- Can be customized for desired usage - whether as a sling, post-operative treatment or for rehabilitation
- Can be worn during competition

Features:

- 1/8" VeLPlush™ neoprene
- Two 1" x 3" hook straps originating from the bicep area to attach to body of brace
- Two additional straps and 1" squares of hook for additional support placement
- Fits right or left shoulder - male or female

• **L3660 approved**
Patent 7,785,281



K195 "THE WINNER™"

Indications:

- Patella tendinitis
- Osgood-Schlatter disease
- Chondromalacia
- Patella-Femoral Syndrome
- Jumper's or runner's knee

Features:

- 1/8" UBL 1-side, non-slip 1-side neoprene
- Moldable insert for desired shape
- Universal up to size 2xlarge



G100 GROIN HAMSTRING STABILIZER

Indications:

- Groin strain
- Hamstring strain
- Quadricep strain

Features:

- Hip Spica design prevents migration
- 2 - 6" Elastic Straps provide superior support
- Soft VeLPlush material
- Specify Left or Right

Product Variation:

G200 Bilateral support



T1 THIGH SUPPORT

Indications:

- Hamstring injuries
- Quadricep pulls

Features:

- 1/8" nylon 4-way stretch neoprene
- 12" sleeve
- Oval pressure pad



Proudly made in the USA

LUMBAR SUPPORTS



L1 ELASTIC DOUBLE-SIDE PULL LUMBAR SUPPORT

Indications:

- Preventative measures
- Lower back pain
- Post-operative care
- Facet Syndrome
- Muscle spasms

Features:

- 9" wide elastic belt
- Two 4" wide adjustable side pulls that overlap in front to increase intra-abdominal support
- 1/8" neoprene pocket for optional insert
- **L0625 Approved**
- Available in black or tan



L9 INSERT

Features:

- 3/16" neoprene
- 1/4" sewn-in pad
- Fits L1 pocket



L10 INSERT

Features:

- Heat moldable
- Fits L1 pocket



L5N ELASTIC DOUBLE-SIDE PULL LUMBAR SUPPORT

Features:

- 12" wide elastic belt
- Two 4" wide adjustable side pulls that overlap in front to increase intra-abdominal support
- 1/8" nylon 4-way stretch neoprene pocket for optional insert
- **L0625 Approved**
- Available in black only



L11 INSERT

Features:

- 11" height
- Heat moldable
- Fits L5N pocket



LC10 KOOLFLEX™ LUMBAR SUPPORT

Features:

- Breathable, hypo-allergenic **Koolflex™** posterior panel
- Latex-Free, Elastic anterior panels with hook and loop closure
- Four removable pre-contoured malleable aluminum posterior stays
- Two reinforced circumferential straps
- Heavy duty reinforced loops
- Anterior hook and loop Velcro
- Available in black only

FRONT



BACK

STAYS & HINGES

Pictured Left to Right

S -Spiral stay w/ dipped caps

PC-Polycentric, stainless steel 1820 hinge w/ hyperextension stop

MP -Multi-Position 1832 hinge w/ 0°, 15°, 30°, 45° extension stops

HT -Polycentric flip-top, high-tech 1832 hinge w/ easy set flexion and extension stops w/ aircraft aluminum arms

NOS - Uncovered 1832 hinge with flexion and extension stops



C-ME™ PROSTHETIC COVER SLEEVES



C-ME™ PROSTHETIC COVER SLEEVES and LAMINATING PANELS FOR AMPUTEES

Features:

- Made in the U.S.A. of the finest quality of Spandex, Lycra and Nylon materials
- Available in a variety of colors and designs to enhance your attitude and personality
- Looks great w/ short pants; helps long pants slide on better
- Change colors and designs as you do your clothing; helps give you a bright outlook every day
- Many patterns available
- BK Sleeve Cover - helps protect the sleeve and socket

Machine wash gentle cycle, cold water. Air-dry only

Laminating Panels measure 20" x 20"



Size	Sleeve Length	Top Circumference	Bottom Circumference
XS (Pedi)	11" - 15"	11" - 14"	8" - 10"
S/M	17" - 21"	15" - 18"	9" - 12"
M/L	17" - 21"	16" - 20"	10" - 14"
XL	19" - 23"	18" - 24"	11" - 15"

PROSTHETIC SUSPENSION SLEEVES



SSL50



SSL51



SSL52



SSL54



SSL55

SUSPENSION SLEEVES

Features:

Airtight seams hold prosthesis secure

SSL50 - 1/8" 1-part suspension sleeve w/ flexion knee, 2-sides nylon

SSL51 - 1/8" 1-part suspension sleeve w/ flexion knee, 1-side
sharkskin neoprene

SSL52 - 1/8" 2-part suspension sleeve w/ flexion knee & top strap

SSL54 - 1/8" 3-part suspension sleeve w/ flexion knee

SSL55 - 1/8" 3-part suspension sleeve w/ lycra middle

Measure circumference at
5" above the knee

Sleeves shown and stocked
in Tan

All suspension sleeves
available tan, golden tan, black
and brown



Proudly made in the USA

OAB PART NUMBERS

Item Number	Description
OABL-XS	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, XS
OABL-S	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, S
OABL-M	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, M
OABL-LG	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, LG
OABL-XL	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, XL
OABL-2XL	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, 2XL
OABL-3XL	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, 3XL
OABL-4XL	"New Option OA" Brace, Hinge Fits Left Medial or Right Lateral OA, 4XL
OABR-XS	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, XS
OABR-S	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, S
OABR-M	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, M
OABR-LG	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, LG
OABR-XL	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, XL
OABR-2XL	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, 2XL
OABR-3XL	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, 3XL
OABR-4XL	"New Option OA" Brace, Hinge Fits Right Medial or Left Lateral OA, 4XL

OAW PART NUMBERS

Item Number	Description
OAWL-XS	"Wrap Around OA" Brace-Hinge Fits Left Knee Medial or Right Knee Lateral OA, XS
OAWL-SM	"Wrap Around OA" Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, S
OAWL-MD	"Wrap Around OA" Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, M
OAWL-LG	"Wrap Around OA" Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, LG
OAWL-XL	"Wrap Around OA" Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, XL
OAWL-2XL	"Wrap Around OA" Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, 2XL
OAWL-3XL	"Wrap Around OA" Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, 3XL
OAWL-4XL	"Wrap Around OA" Brace, Hinge Fits Left Knee Medial or Right Knee Lateral OA, 4XL
OAWR-XS	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, XS
OAWR-SM	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, S
OAWR-MD	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, M
OAWR-LG	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, LG
OAWR-XL	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, XL
OAWR-2XL	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, 2XL
OAWR-3XL	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, 3XL
OAWR-4XL	"Wrap Around OA" Brace, Hinge Fits Right Knee Medial or Left Knee Lateral OA, 4XL



FITTING GUIDE

NOTE: There will be an up-charge for all custom-made products.

XS - 3XL Professional Price

4XL \$10.00 up charge fee

KNEE SUPPORTS

Measure around the center of the knee with leg extended.

	Mid-patella	6" above knee
XSMALL	12" - 13"	16" - 19"
SMALL	13" - 14"	17" - 20"
MEDIUM	14" - 15"	18" - 21"
LARGE	15" - 16"	19" - 22"
XLARGE	16" - 18"	20" - 24"
2XLARGE	18" - 20"	22" - 26"
3XLARGE	20" - 22"	24" - 28"
4XLARGE	22" - 24"	26" - 30"
CUSTOM	Please provide measurements	

WRIST SUPPORTS

Measure around wrist

XSMALL	4" - 5"
SMALL	5" - 6"
MEDIUM	6" - 7"
LARGE	7" - 8"
XLARGE	8" - 9"
2XLARGE	9" - 10"
3XLARGE	10" - 11"
CUSTOM	11" + UP

ELBOW SUPPORTS

Measure the largest circumference of the forearm.

X-SMALL	8" - 9"
SMALL	9" - 10"
MEDIUM	10" - 11"
LARGE	11" - 12"
X-LARGE	12" - 14"
XX-LARGE	14" - 16"
XXX-LARGE	16" - 18"
CUSTOM	18" + UP

A50 ANKLE SUPPORT

By shoe size

Size may vary according to ankle circumference.

	Men's	Women's
XSMALL	5 - 6	7 - 8
SMALL	6 - 7	8 - 9
MEDIUM	7 - 9	9 - 11
LARGE	10 - 12	12 - 14
XLARGE	12 - 14	14 - 16
2XLARGE	14 - 16	

A30-WOOTEN UHL

By men's shoe size.

X-SMALL	4 - 6
SM/MED	6.5 - 11
LG/XL	11.5 - 13.5
XXL	14 + UP

S5 SHOULDER

Chest circumference Bicep circumference

XSMALL	22" - 25"	10" - 13"
SMALL	25" - 28"	11" - 14"
MEDIUM	28" - 31"	11.5" - 14.5"
LARGE	31" - 35"	12" - 14.75"
XLARGE	36" - 40"	12.5" - 15.5"
2XLARGE	38" - 43"	13" - 16"
3XLARGE	42" - 48"	13.5" - 16.5"
4XLARGE	47" - 52"	14.5" - 17.5"
CUSTOM	Please provide measurements	

CALF SUPPORT

Measure largest circumference of the lower leg.

X-SMALL	12" - 13"
SMALL	13" - 14"
MEDIUM	14" - 15"
LARGE	15" - 16"
X-LARGE	16" - 17"
XX-LARGE	17" - 18"
XXX-LARGE	18" - 19"
CUSTOM	19" + UP

G100 GROIN HAMSTRING STABILIZER

Measure waist circumference

XS	28"-32"
SM/MD	32"-36"
LG/XL	36"-40"
XXL	40"-44"
CUSTOM	44" & up

THIGH SUPPORT

Measure mid-thigh circumference.

X-SMALL	18" - 19"
SMALL	19" - 20"
MEDIUM	20" - 22"
LARGE	22" - 25"
X-LARGE	25" - 28"
XX-LARGE	28" - 30"
XXX-LARGE	31" - 34"
CUSTOM	34" + UP

LUMBAR SUPPORT

Measure around waist.

X-SMALL	22" - 26"
SMALL	26" - 30"
MEDIUM	30" - 34"
LARGE	34" - 38"
X-LARGE	38" - 42"
XX-LARGE	42" - 46"
XXX-LARGE	46" - 50"
XXXX-LARGE	50" - 54"
CUSTOM	54" + UP

LC10 Special Sizing Guide

Measure around waist.

XSMALL	25" - 32"
SMALL	28" - 35"
MEDIUM	32" - 41"
LARGE	34" - 44"
XLARGE	40" - 48"
2XLARGE	42" - 54"
3XLARGE	44" - 59"

SUSPENSION SLEEVES

Measure circumference at 5" above the knee

Size 000	8" - 9.5"
Size 00	9.5" - 10.5"
Size 0	10.5" - 12.5"
Size 1	11.5" - 13"
Size 1.5	12.5" - 14"
Size 2	14" - 16"
Size 2.5	15" - 17"
Size 3	16" - 18"
Size 3.5	17" - 18.5"
Size 4	18.5" - 20"
Size 5	21" - 23"



Proudly made in the USA



CUSTOMS and PRIVATE LOGO

All products are in stock and shipped no later than the following business day, with the exception of very large orders.

NOS has a Private Label Logo program available. Check with Customer Service and Sales for details.

Exclusive CUSTOMS Program - see prices and details below:

Knee and Elbow Sleeves, Wrists and Suspension Sleeves - Custom Fee \$35.00

Hinged Knee & Elbows, Lumbars, Groin and Shoulder, and any products not mentioned above. - Custom Fee \$75.00

All made to order Custom Braces are **NON-RETURNABLE**. Shipment is normally within five (5) working days. Prices for Custom Braces are priced as listed in the current New Options Sports Professional Price List plus the non-discountable custom charge.

Please refer to New Options Sports' Professional Price List for custom measuring guides.

Color change available for some materials and braces. All color change orders are **NON-RETURNABLE**. Shipment is normally within five (5) working days.

New Options Sports, 1850 Diplomat Drive #100, Farmers Branch, TX 75234
PH: 800.872.5488 or 214.638.6422 FX: 800.455.5488 or 214.638.6425
Questions: info@newoptionssports.com Orders: orders@newoptionssports.com