

## TECHFORM® CASTING TAPE

Techform Casting Tape was designed by and for orthopaedic casting professionals who look for a superior handling experience and a smooth, strong finished cast. It incorporates a patented rub-responsive formula that works with our unique fabric to produce the best finished cast you'll ever see.



	Techform Premium
Fabric	Fiberglass and polypropylene blend
Resin	Rub-responsive formula
Precision-welded ends	Yes
Fray-resistant cutting	No
Fabric Weight	Light
Finished cast surface	Smooth
Fluffy rolls for easier application	Yes
Core with holes for consistent set time	Yes
Low-profile packaging with finger holes	Yes

### TECHFORM PREMIUM ORDERING INFORMATION (BY COLOR)

WIDTH	WH	BK	OR	RD	PK	HP	LB	DB	PR	GN	RB
1.5"	154WH										
2"	204WH	204BK	204OR	204RD	204PK	204HP	204LB	204DB	204PR	204GN	204RB
3"	304WH	304BK	304OR	304RD	304PK	304HP	304LB	304DB	304PR	304GN	304RB
4"	404WH	404BK	404OR	404RD	404PK	404HP	404LB	404DB	404PR	404GN	404RB
5"	504WH										

White (WH)	Black (BK)	Orange (OR)	Red (RD)	Pink (PK)
Hot Pink (HP)	Li-Blue (LB)	Dk Blue (DB)	Purple (PR)	Green (GN)

- 10 individually-pouched, 4-yard rolls per box
- Rainbow box (RB) contains 2 rolls each of 5 different colors (light/dark blue, hot pink, red & black)

