

ÖSSUR®

FORM FIT®

THE CHOICE OF EXPERTS

GSA / FSS CONTRACT NO:
V797D60740



Form Fit Walker

Sometimes, our bodies need assistance. For decades, specialists have trusted Form Fit supports, a complete range of quality products by Össur. A global leader in innovative orthopaedic products, Össur's biomechanical expertise has helped millions live Life Without Limitations.

That's why Form Fit is the Choice of Experts.



FORMFIT®

THE CHOICE OF EXPERTS

FORM FIT WALKER

The Form Fit® Walker Air provides immobilization for patients suffering from soft tissue injury; grades 2 & 3 strains and sprains, and stable fractures.

FORM FIT WALKER, AIR

The patented pneumatic system is designed to increase stability while decreasing pain and swelling. The intuitive pump and clearly marked decompression valve make air/pressure adjustments easy.

Features and Benefits

- High and low top versions available
- Integrated pneumatic system (Air version only)
- Low, wide rocker bottom promotes a stable, natural gait
- Shock absorbing sole for reduced impact at heel strike
- Breathable, seamless liner

Indications

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative, trauma or rehab use



GSA / FSS CONTRACT NO: V797D60740

FORM FIT WALKER

Size	Men's Shoe	Women's Shoe	Standard	Low Top
XSmall	2-4	3.5-5.5	W0200BLK	W0300BLK
Small	4.5-7	6-8	W0400BLK	W0500BLK
Medium	7.5-10.5	8.5-11.5	W0600BLK	W0700BLK
Large	10.5-12.5	11.5-13.5	W0800BLK	W0900BLK
XLarge	12.5+	13.5+	W1000BLK	W1100BLK

Replacement Softgoods

Size	Standard	Low Top
XSmall	WKSG200BLK	WKSG0300BLK
Small	WKSG400BLK	WKSG0500BLK
Medium	WKSG600BLK	WKSG0700BLK
Large	WKSG800BLK	WKSG0900BLK
XLarge	WKSG1000BLK	WKSG1100BLK

FORM FIT WALKER, AIR

Size	Men's Shoe	Women's Shoe	Standard	Low Top
XSmall	2-4	3.5-5.5	A-W0200BLK	A-W0300BLK
Small	4.5-7	6-8	A-W0400BLK	A-W0500BLK
Medium	7.5-10.5	8.5-11.5	A-W0600BLK	A-W0700BLK
Large	10.5-12.5	11.5-13.5	A-W0800BLK	A-W0900BLK
XLarge	12.5+	13.5+	A-W1000BLK	A-W1100BLK

Replacement Softgoods

Size	Standard	Low Top
XSmall	12349	12356
Small	10675	12357
Medium	10678	12358
Large	10681	12359
XLarge	12347	12360

ORDERS: fax 720-221-0677 | orders@frontlinemedicalmil.com | phone 866-681-0444



Frontline Medical Services
Toll free: 866-681-0444
www.frontlinemedicalmil.com



Serve proud. Live well.™

GSA / FSS CONTRACT NO: V797D60740