

## UNLOADER ONE<sup>®</sup>

The Unloader One is a comfortable, lightweight, low-profile brace that provides excellent suspension and pain relief. This lightweight brace preserves the clinically proven Össur Unloader 3-Points of Leverage design by incorporating Dual Dynamic Force Straps that are easily fine-tuned with the SmartDosing<sup>®</sup> dial. Color-coded Quick Fit Buckles simplify application and removal, while flexible shells and breathable tibial liners, coated with Össur Sensil<sup>®</sup> Silicone, improve comfort and reduce migration, helping to ensure patient compliance. Recent research demonstrates that the Unloader One reduces pain, improves function and decreases the use of pain medication.



### INDICATIONS FOR USE

- Moderate to severe unicompartmental osteoarthritis

### FEATURES AND BENEFITS

- Dual Dynamic Force Straps designed to reduce pressure on sensitive bony surfaces of the knee, resulting in reduction of pain
- Tool-free SmartDosing dial enables fine-tuned pain management and can be adjusted even when wearing the brace under clothing
- Color coded quick fit buckles enable easy application and removal
- Proprietary Sensil<sup>®</sup> silicone liners designed to reduce migration
- Flexible, low-profile upright contours to patient anatomy and allows use underneath clothing
- Quick Connect Strapping for tool-free removal and replacement
- Lightweight frame is perforated for flexibility and ventilation with custom paint options for the Custom version
- Unloader One is available in pre-fabricated and Custom version

### REQUIRED MEASUREMENTS

- Circumference measurement: 6" (15 cm) below mid-patella



## UNLOADER ONE®

### ORDERING INFORMATION

#### UNLOADER ONE WITH SMARTDOSING

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2 – 14" (31 – 36 cm)	B-240618712	B-240518712	B-240628712	B-240528712
Medium	14 – 16.75" (36 – 43 cm)	B-240618713	B-240518713	B-240628713	B-240528713
Large	16.75 – 20" (43 – 51 cm)	B-240618714	B-240518714	B-240628714	B-240528714
X-Large	20 – 24" (51 – 61 cm)	B-240618715	B-240518715	B-240628715	B-240528715

#### UNLOADER ONE SHORT WITH SMARTDOSING

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	9.8 – 12.2" (25 – 31 cm)	B-241618711	B-241518711	B-241628711	B-241528711
Small	12.2 – 14" (31 – 36 cm)	B-241618712	B-241518712	B-241628712	B-241528712
Medium	14 – 16.75" (36 – 43 cm)	B-241618713	B-241518713	B-241628713	B-241528713
Large	16.75 – 20" (43 – 51 cm)	B-241618714	B-241518714	B-241628714	B-241528714
X-Large	20 – 24" (51 – 61 cm)	B-241618715	B-241518715	B-241628715	B-241528715

#### UNLOADER ONE WITH RATCHET

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2 – 14" (31 – 36 cm)	B-240619712	B-240519712	B-240629712	B-240529712
Medium	14 – 16.75" (36 – 43 cm)	B-240619713	B-240519713	B-240629713	B-240529713
Large	16.75 – 20" (43 – 51 cm)	B-240619714	B-240519714	B-240629714	B-240529714
X-Large	20 – 24" (51 – 61 cm)	B-240619715	B-240519715	B-240629715	B-240529715

#### UNLOADER ONE SHORT WITH RATCHET

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	9.8 – 12.2" (25 – 31 cm)	B-241619711	B-241519711	B-241629711	B-241529711
Small	12.2 – 14" (31 – 36 cm)	B-241619712	B-241519712	B-241629712	B-241529712
Medium	14 – 16.75" (36 – 43 cm)	B-241619713	B-241519713	B-241629713	B-241529713
Large	16.75 – 20" (43 – 51 cm)	B-241619714	B-241519714	B-241629714	B-241529714
X-Large	20 – 24" (51 – 61 cm)	B-241619715	B-241519715	B-241629715	B-241529715

Standard version is recommended for patients between 5'5" and 6'4" (1.65 and 1.92 meters tall).

Short version is recommended for patients under 5'5" (1.65 meters tall).

Hinge is placed on the side of the affected compartment

GSA / FSS CONTRACT NO: V797D60740

ORDERS: fax 720-221-0677 | orders@frontlinemedicalmil.com | phone 866-681-0444



Frontline Medical Services  
Toll free: 866-681-0444  
www.frontlinemedicalmil.com

GSA / FSS CONTRACT NO: V797D60740



Serve proud. Live well.™