

## Freestyle™ OA Knee Brace

The Freestyle OA provides off-loading of the medial or lateral compartment for patients with mild to moderate unicompartmental osteoarthritis. The simple wraparound design allows patients to apply and remove the brace with ease. Comfortable, breathable material provides focused compression around the knee joint. The single-sided low profile hinge encourages patient compliance, while the innovative thumbwheel dial allows precision unloading adjustments without the use of a tool.

### Common Examples of Use

- Mild to unicompartmental osteoarthritis

### Features

- Lightweight, low profile brace fits under most clothing

- Soft, comfortable material promotes patient compliance
- Simple, intuitive wraparound design provides compression
- Adjustable off-loading – no tools required
- Finger pockets aid in application
- Medial or lateral options available



Description	Part #
Freestyle OA Medial, Left	1171X
Freestyle OA Medial, Right	1172X
Freestyle OA Lateral, Left	1176X
Freestyle OA Lateral, Right	1177X

X = See size chart.

CE Marking

\*When ordering extension kits for the Freestyle OA Lateral, you will need to order the kit for the opposite leg so the application number sequence matches the brace.

Description	Part #
<b>Accessories</b>	
Freestyle OA Condyle Pad	100338-000
Freestyle OA Thigh Extension Kit, Left*	117000-100
Freestyle OA Thigh Extension Kit, Right*	117000-200
Freestyle OA Calf Extension Kit, Left*	117001-100
Freestyle OA Calf Extension Kit, Right*	117001-200

Freestyle OA Size Chart			
Size #	Size	6" Above Mid Patella	6" Below Mid Patella
X = 2	S	15" - 18" (38-46 cm)	12" - 14" (30-36 cm)
X = 3	M	18" - 21" (46-53 cm)	14" - 16" (36-41 cm)
X = 4	L	21" - 24" (53-61 cm)	16" - 18" (41-46 cm)
X = 5	XL	24" - 27" (61-69 cm)	18" - 20" (46-51 cm)
X = 6	XXL	27" - 30" (69-76 cm)	20" - 23" (51-58 cm)

Circumference measured 6" (15 cm) above mid patella and 6" (15 cm) below mid patella.

ORDERS: fax 720-221-0677 | orders@frontlinemedicalmil.com | phone 866-681-0444



Frontline Medical Services  
Toll free: 866-681-0444  
www.frontlinemedicalmil.com



Serve proud. Live well.™

