

## Legacy Thruster

The Legacy Thruster is a small, lightweight single-upright brace with the hinge on the affected side of the knee. The straps are pulled to open the joint space and provide relief. This design addresses OA conditions only, not associated with ligament instability.

### Common Examples of Use

- Unicompartmental OA
- Meniscal instabilities

### Features

- Low-profile frame
- Dial allows incremental force adjustments
- Hinge is on affected side of knee, pulling vs. pushing to help relieve the symptoms of OA
- Unloads in a 4-point load system
- Distal calf section flexes with the muscles of the leg
- Up to 22° of correction can be applied
- For moderate to severe OA patients who require the brace for all-day use



Custom color shown

Description	Part #
Legacy Thruster, Custom	TA409000
Legacy Thruster, Custom, Short	TA419000
Legacy Thruster, Lateral, XS - XXL	TA408YXX
Legacy Thruster, Medial, XS - XXL	TA409YXX
Legacy Thruster, Lateral Short, XS - XXL	TA418YXX
Legacy Thruster, Medial Short, XS - XXL	TA419YXX
Legacy Thruster, Lateral Athletic, XS - XXL	TA428YXX
Legacy Thruster, Medial Athletic, XS - XXL	TA429YXX
Legacy Thruster, Lateral Short Athletic, XS - XXL	TA438YXX
Legacy Thruster, Medial Short Athletic, XS - XXL	TA439YXX

XX = See size chart. Y = Left 1, Right 2  
CE Marking

Description	Part #
<b>Soft Goods Kit</b>	
Custom	TA460000
Medial Left / Lateral Right, XS - XXL	TA4601XX
Medial Right / Lateral Left, XS - XXL	TA4602XX
Short, Medial Left / Lateral Right, XS - XXL	TA4611XX
Short, Medial Right / Lateral Left, XS - XXL	TA4612XX
Athletic, Medial Left / Lateral Right, XS - XXL	TA4621XX
Athletic, Medial Right / Lateral Left, XS - XXL	TA4622XX
Athletic Short, Medial Left / Lateral Right, XS - XXL	TA4631XX
Athletic Short, Medial Right / Lateral Left, XS - XXL	TA4632XX

The same kit works for a medial left and a lateral right brace.

The same kit works for a medial right and a lateral left brace.

Description	Part #
<b>Accessories</b>	
Calf Pad Kit ML/LR <sup>1,2</sup>	TA4651XX
Calf Pad Kit MR/LL <sup>1,2</sup>	TA4652XX
Calf Pad Kit ML/LR, Short <sup>1,2</sup>	TA4671XX
Calf Pad Kit MR/LL, Short <sup>1,2</sup>	TA4672XX
Calf Strap Kit <sup>1,2</sup>	TA4610XX
Thigh Pad Kit ML/LR <sup>2</sup>	TA4641XX
Thigh Pad Kit MR/LL <sup>2</sup>	TA4642XX
Thigh Pad Kit ML/LR, Short <sup>2</sup>	TA4661XX
Thigh Pad Kit MR/LL, Short <sup>2</sup>	TA4662XX
Thigh Strap Kit <sup>2</sup>	TA4600XX

<sup>1</sup>Also available in XXS (XX = 00).

<sup>2</sup>Not available in 3XL.

### Legacy Thruster Size Chart

Size #	Size	Thigh Circumference	Standard Calf Circ.	Athletic Calf Circ.
XX = 01	XS	15" - 17.5"	12.5" - 14"	10.5" - 12"
XX = 03	S	17.5" - 20"	14" - 15.5"	12" - 13.5"
XX = 05	M	20" - 22.5"	15.5" - 17"	13.5" - 15"
XX = 07	L	22.5" - 25"	17" - 18.5"	15" - 16.5"
XX = 09	XL	25" - 27.5"	18.5" - 20"	16.5" - 18"
XX = 11	XXL	27.5" - 29"	20" - 21.5"	18" - 19.5"

Thigh circumference measure 6" (15 cm) above mid-patella. Calf circumference measure 6" (15 cm) below mid-patella.

ORDERS: fax 720-221-0677 | orders@frontlinemedicalmil.com | phone 866-681-0444



Frontline Medical Services  
Toll free: 866-681-0444  
www.frontlinemedicalmil.com



Serve proud. Live well.™