

## Polar Care Cube™

The Polar Care Cube Cold Therapy System is simple and reliable. Its easy to use, compact design makes it great for clinic, hospital, and home use. With the addition of ice and water, your patients will enjoy 6-8 hours of effortless cold therapy. Proper use requires an insulation barrier between the pad and the patient's skin. Federal law restricts this device to sale by or on the order of a medical professional.

### Common Examples of Use

- Post-operative
- Arthroscopic procedures
- Reconstructive procedures
- Plastic surgery
- General surgery
- Post trauma
- Chronic pain
- Physical therapy

### Features

- Physicians and patients benefit from reduced pain and swelling—follow-up visits focus on recovery, rather than pain management
- Whisper-quiet pump technology allows patients to rest and recover in comfort
- Durable, powerful construction provides peace of mind
- Broad selection of premium WrapOn Pads conform to the affected area
- Sterile Polar Dressings available for each part



Description	Part #
Polar Care Cube	10701
<b>Cube Combo Units with WrapOn Pad (includes Cube and 1 single-patient use pad)</b>	
PC Cube with WrapOn Knee Pad	10705
Sterile Polar Dressing, Knee L	04708
PC Cube with WrapOn L Knee Pad	10706
Sterile Polar Dressing, Knee L	04708
PC Cube with WrapOn XL Knee Pad	10707
Sterile Polar Dressing, Shoulder	04908
PC Cube with WrapOn Ankle Pad	10708
Sterile Polar Dressing, Back, Ankle, Wrist	09810
PC Cube with WrapOn Multi-Use XL Pad with Long Stem (16" Long)	10709
Sterile Polar Dressing, XL Multi-Use	02344
PC Cube with WrapOn XL Multi-Use Pad	10710
Sterile Polar Dressing, XL Multi-Use	02344
PC Cube with WrapOn Shoulder Pad	10711
Sterile Polar Dressing, Shoulder	04908
PC Cube with WrapOn XL Shoulder Pad	10712
Sterile Polar Dressing, XL Shoulder	04918
PC Cube with WrapOn Hip Pad	10713
Sterile Intelli-Flo Dressing, Hip	10680

Note: Matching Sterile Polar Dressing listed under each pad. Body part indicated for Sterile Polar Dressing may or may not correlate with body part in pad name, but the sizes have been paired appropriately.

Description	Part #
<b>Bulk Packs</b>	
Polar Care Cube, 4-Pack	28600
PC Cube Combo, WrapOn XL Multi-Use, 4-Pack	28601
PC Cube Combo, WrapOn Knee, 4-Pack	28602
PC Cube Combo, WrapOn Shoulder, 4-Pack	28603
<b>PC Cube Combo Units with Polar Pad (includes Cube and 1 single-patient use pad)</b>	
PC Cube with Polar Knee / Shoulder Pad	10702
Sterile Polar Dressing, Knee / Shoulder	02328
PC Cube with Polar Multi-Use Pad	10704
Sterile Polar Dressing, Multi-Use	02348
<b>Accessories</b>	
Replacement Power Supply	10698
Compression Strap (60" long)	10670
PC Y Adaptor	07640
Polar Care Bed Hanger	000038-000
4ft Extension Tube	09915

ORDERS: fax 720-221-0677 | orders@frontlinemedicalmil.com | phone 866-681-0444



Frontline Medical Services  
Toll free: 866-681-0444  
www.frontlinemedicalmil.com



Serve proud. Live well.™