

JEANIE RUB® VARIABLE SPEED MASSAGER

- Deep, penetrating massage with a twist of the switch
- Helps reduce pain, muscle tension, and stiffness
- Variable speeds provide soothing massage or more invigorating massage
- An industry standard among professional massage therapists
- Helps stimulate the circulatory system
- Fosters healing of strained muscles and sprained ligaments



Available from:



www.frontlinemedicalmil.com



THERAPEUTIC RELIEF

The Jeanie Rub® Variable Speed Massager delivers 1100 - 4600 rpm for a deep, penetrating massage. Help reduce pain, muscle tension, and stiffness at the twist of the switch. The orbital action of the massager provides a softer and deeper massage than "percussion" type massagers. Lower speeds deliver a soothing massage while faster speeds are more invigorating.

HIGH QUALITY CONSTRUCTION

This product is built with a cast aluminum construction and has a red and black powder coat finish that looks great and resists abrasions, ensuring years of use. Jeanie Rub massagers are built in the USA to ANSI Standards and are CSA certified, ensuring they have been evaluated to meet or exceed consumer safety standards. Pad measures 4" x 9" (10cm x 23cm). Jeanie Rub accessories sold separately. 12' grounded power cord, 120V AC Nominal power source required.

SPECIFICATIONS

Product Selection:
PRO-3401 Jeanie Rub®
Variable Speed Massager

Product Weight: 7 lbs
Dimensions: 10.25" x 8.25"

Case Dimensions:
19.125" x 10.875" x 8.5"
Case Weight: 28 lbs
Case Quantity: 4

Latex: This product is not made with natural rubber latex.

JEANIE RUB® VARIABLE SPEED MASSAGER ACCESSORIES

Professional package includes:

- Jeanie Rub Variable Speed Massager
- Paraspinal Accessory
- Extremity Accessory
- Carrying Bag
- Fleece Cover

*Professional Package
PRO-3405*



*Paraspinal Accessory
ACC-896*

*Extremity Accessory
ACC-897*



*Wall Hanger
ACC-884*



*Fleece Cover
ACC-882*



*Sheepskin Cover
ACC-886*



*Carrying Bag
ACC-885*



*Terry Cover
ACC-893*



*Sure-Chek Cover
ACC-892*

SPECIFICATIONS

Product Selection:	Dimensions:	Wt:	Case Qty:	Case Wt:
PRO-3405: Professional Package	13.5" x 10.25" x 12.5"	11 #	1	11 #
ACC-896 Paraspinal Accessory	10.25" x 5" x 2"	.15 #	3	2 #
ACC-897 Extremity Accessory	10.25" x 5" x 4.5"	.2 #	6	2 #
ACC-884 Wall Hanger	14" x 5" x 5"	1.5 #	1	3 #
ACC-882 Fleece Cover	9.5" x 5.5" x .5"	.1 #	12	4 #
ACC-886 Sheepskin Cover	10" x 8" x 1"	.3 #	6	3 #
ACC-885 Carrying Bag	12" x 8" x 5"	1 #	6	2 #
ACC-893 Terry Cover	10" x 7"	.2 #	3	2 #
ACC-892 Sure-Chek Cover	9" x 4"	.15 #	3	2 #

Available from:



Serve proud. Live well.™

www.frontlinemedicalmil.com