

SWEDE-O® INNER LOK 8®

- Easy to apply and adjust using internal figure eight straps
- Simple removal with quick release strap
- Stays tight and secure with exclusive lacing system and top strap
- Durable, lightweight ballistic nylon is ideal for sports
- Full elastic back reduces irritation and potential for blistering
- Bilateral design fits right or left ankle
- Industry leading 1 year warranty, satisfaction guarantee



Available from:



www.frontlinemedicalmil.com



THE BRACE THAT GIVES YOU INNER STRENGTH

The Swede-O® Inner Lok 8® features a patented internal figure eight strapping system which assures proper strap placement and simulates professional taping.

The quick release strap allows you to easily remove the brace by loosening the internal strapping system. The exclusive dual eyelet lacing system is offset to create a "pulley effect", keeping the brace tighter longer. The top strap covers the laces and locks the figure eight straps in position.

The full elastic back is comfortable against the Achilles tendon and helps reduce irritation and potential for blistering. The breathable tongue helps keep the foot cool and comfortable. The Inner Lok 8 is backed by an industry leading 1 year warranty and constructed from a durable ballistic nylon which is comfortable thin, and lightweight.

INDICATIONS FOR USE

Use the Inner Lok 8 after injury or to protect your ankle from future injury. This ankle brace helps prevent and relieve aches, pains, sprains, and other ankle conditions.

SPECIFICATIONS

Product Selection:
AKL-6335 Swede-O® Inner Lok 8®

Product Weight: .35 lb
Packaging Dimensions: 13" x 9" x 1"

Case Weight: 6 lbs

Case Dimensions: 11" x 14" x 13"
Case Quantity: 12

Materials: Polyester, nylon, polyurethane foam, brass

Patent: 6,398,750 & 6,652,474

Latex: This product is not made with natural rubber latex.

Sized by US shoe size

MODEL #	SIZE	MEN'S	WOMEN'S
AKL-6335-BK-1XS	XS	3-5 1/2	4-6 1/2
AKL-6335-BK-SML	SM	6-7 1/2	7-9 1/2
AKL-6335-BK-MED	MED	8-10 1/2	10-11 1/2
AKL-6335-BK-LRG	LRG	11-12 1/2	12-13 1/2
AKL-6335-BK-1XL	XL	13-14 1/2	14-15 1/2