

## THERMAL VENT® PLANTAR DR™

- Nighttime plantar fasciitis relief
- Helps relieve heel pain
- Warmth, gentle stretching, and light compression promotes healing
- Lightweight and portable
- Bilateral design fits both right and left foot



### THERMAL TECHNOLOGY

Engineered with MVT<sup>2</sup>:  
MicroVentilated Thermal Technology™

Thermal Vent Supports are a valuable aid to help rehabilitate injured tissue. These products provide warmth and heat therapy to help increase pliability and elasticity of muscles and tendons, making them more responsive to therapy and exercise.



Available from:



Serve proud. Live well.™

www.frontlinemedicalmil.com



### COMFORTABLE PLANTAR FASCIITIS SUPPORT TO RELIEVE HEEL PAIN

The Swede-O® Thermal Vent® Plantar DR™ is designed to help relieve pain and discomfort from plantar fasciitis, Achilles tendinitis, and other conditions.

The non-elastic strap pulls the ankle into a slight dorsiflexion position and extends the toes to gently stretch the plantar fascia; passively stretching the calf and Achilles tendon while resting or sleeping to help alleviate morning pain and discomfort. The gentle stretch helps reduce muscle contracture, inflammation, and associated pain. Non-slip safety sole helps prevent slipping.

### SPECIFICATIONS

#### Product Selection:

BRE-6340 Swede-O® Thermal Vent® Plantar DR™

**Product Weight:** .3 lb

**Dimensions:** 13" x 9" x 1.5"

**Case Weight:** 5 lbs

**Case Dimensions:** 11" x 14" x 13"

**Case Quantity:** 12

**Materials:** Outer lining: Nylon, Mid Layer: Neoprene, Inner Lining: Nylon.

**Latex:** This product is not made with natural rubber latex.

Sized by US shoe size

MODEL #	SIZE	MEN'S	WOMEN'S
BRE-6340-GR-1XS	XS	3-5	4-6
BRE-6340-GR-SML	SM	5 1/2-7	6 1/2-9
BRE-6340-GR-MED	MED	7 1/2-10	9 1/2-11
BRE-6340-GR-LRG	LRG	10 1/2-12	11 1/2-13
BRE-6340-GR-1XL	XL	12 1/2-14	13 1/2-15

Your Thermal Vent Support should be a snug and comfortable fit.  
Not too tight, not too loose. If in doubt, try next larger size.