

# AIR WALKER

Orders may be placed through the DLC Remote Order Entry System (ROES) online ordering application. For assistance, contact DLC Customer Service Support at 303-273-6200.

Contract No: 36C791-20-D-0009



The Össur Logo, trademarks and images are used with permission of Össur hf.



Frontline Medical Services  
Toll Free: 866-681-0444  
[www.frontlinemedicalmil.com](http://www.frontlinemedicalmil.com)



Serve proud. Live well.™



# AIR WALKER

The Air Walker provides immobilization for patients suffering from soft tissue injury; grades 2 & 3 strains and sprains, and stable fractures. The patented pneumatic system is designed to increase stability while decreasing pain and swelling. The intuitive pump and clearly marked decompression valve make air/pressure adjustments easy.

## FEATURES AND BENEFITS

- High and low top versions available
- Integrated pneumatic system
- Low, wide rocker bottom promotes a stable, natural gait
- Shock absorbing sole for reduced impact at heel strike
- Breathable, seamless liner

## INDICATIONS

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative, trauma or rehab use



## ORDERING INFORMATION

Part Number	Description
AW0600BLKC	Centec Air Walker, Medium
AW0800BLKC	Centec Air Walker, Large
AW1000BLKC	Centec Air Walker, XLarge

Part Number	Description
AW0700BLKC	Centec Low Top Air Walker, Medium
AW0900BLKC	Centec Low Top Air Walker, Large
AW1100BLKC	Centec Low Top Air Walker, XLarge



Frontline Medical Services  
Toll Free: 866-681-0444  
www.frontlinemedicalmil.com



Serve proud. Live well.™



PROUD DISTRIBUTOR OF