

COMPRESSOGRIP | PROSTHETIC SHRINKER

Available through VA Denver Logistics Center (DLC)

Contract No: 36C79120D0009

Orders may be placed through the DLC Remote Order Entry System (ROES) online ordering application. For assistance, contact DLC Customer Service Support at 303-273-6200.

Compressogrip B/K Elastic Shrinker - Transtibial

The original tubular compression prosthetic shrinker used for control or prevention of edema following amputation surgery, for limb volumetric maintenance, to minimize hypertrophic scarring and to give compressive support for pain relief.

FEATURES

- Ease of application by patient
- Clinically tested - wearer accepted
- Easily adjusted to varying lengths
- Assures controlled compression
- Available in White, Beige and Grey X-STATIC®
- Grey X-STATIC® has all the benefits of the original Compressogrip Shrinker
- X-STATIC® silver fiber adds permanent antimicrobial benefits to inhibit odor in the garment and effectively wick moisture away from the skin
- Latex-free material
- Machine washable, tumble dry warm

INDICATIONS

- Pre-prosthetic limb shaping
- Volume management of transtibial residual limb

CONTRAINDICATIONS

If rash, blisters, or skin irritations occur, consult your physician and discontinue use.



KNIT-RITE
INCORPORATED



Frontline Medical Services
Toll Free: 866-681-0444
www.frontlinemedicalmil.com



COMPRESSOGRIP | PROSTHETIC SHRINKER

COMPRESSOGRIP B/K ELASTIC SHRINKER SIZING GUIDE

White	Beige	Grey/Grey X-Static	Size	Length	Fits Circumference
2H0S0218	2H0B0218	2H0X0218	#2	18" (45cm)	5"-7" (13-18cm)
2H0S0224	2H0B0224	---	#2	24" (60cm)	5"-7" (13-18cm)
2H0S0518	2H0B0518	2H0X0518	#2.5	18" (45cm)	7"-10" (18-25cm)
2H0S0524	2H0B0524	---	#2.5	24" (60cm)	7"-10" (18-25cm)
2H0S0318	2H0B0318	2H0X0318	#3	18" (45cm)	9"-12" (23-30cm)
2H0S0324	2H0B0324	2H0X0324	#3	24" (60cm)	9"-12" (23-30cm)
2H0S0330	2H0B0330	---	#3	30" (75cm)	9"-12" (23-30cm)
2H0SS318	2H0BS318	---	#3.5	18" (45cm)	10"-13" (25-33cm)
2H0SS324	2H0BS324	---	#3.5	24" (60cm)	10"-13" (25-33cm)
---	---	2H0XS330	#3.5	30" (75cm)	10"-13" (25-33cm)
2H0S0418	2H0B0418	2H0X0418	#4	18" (45cm)	12"-16" (30-40cm)
2H0S0424	2H0B0424	2H0X0424	#4	24" (60cm)	12"-16" (30-40cm)
2H0S0430	2H0B0430	2H0X0430	#4	30" (75cm)	12"-16" (30-40cm)
2H0S0518	2H0B0518	2H0X0518	#5	18" (45cm)	15"-20" (38-50cm)
2H0S0524	2H0B0524	2H0X0524	#5	24" (60cm)	15"-20" (38-50cm)
2H0S0530	2H0B0530	2H0X0530	#5	30" (75cm)	15"-20" (38-50cm)
2H0S0618	2H0B0618	---	#6	18" (45cm)	18"-24" (45-60cm)
2H0S0624	2H0B0624	2H0X0624	#6	24" (60cm)	18"-24" (45-60cm)
2H0S0630	2H0B0630	2H0X0630	#6	30" (75cm)	18"-24" (45-60cm)



HCPC: L8420

Fits Circumference = Measurement of B/K Circumference 2" Below Medial Tibia Tubercle

NOTE: Compressogrip shrinker fabric gives up to 20-25 mmHg compression at 50% stretch. At 100% stretch, the compression is almost doubled. Double layers will almost double compression at both 50% and 100% stretch. Application of more than two layers may exceed the limits of desired therapeutic pressure.



Place orders directly through the VA Denver Logistics Center (DLC). For assistance, contact DLC Customer Service Support at 303-273-6200.



Frontline Medical Services
Toll Free: 866-681-0444
www.frontlinemedicalmil.com



Serve proud. Live well.™