

# MIAMI J<sup>®</sup> SELECT

## Comfort to compliance

Orders may be placed through the DLC Remote Order Entry System (ROES) online ordering application. For assistance, contact DLC Customer Service Support at 303-273-6200.

Contract No: 36C791-20-D-0009



The Össur Logo, trademarks and images are used with permission of Össur hf.



Frontline Medical Services  
Toll Free: 866-681-0444  
[www.frontlinemedicalmil.com](http://www.frontlinemedicalmil.com)



Serve proud. Live well.™



# MIAMI J® SELECT

The Miami J Select cervical collar by Össur offers the widest range of height adjustability in the market with a single, universal design that replicates the comfort, immobilization and quality you expect from a Miami J product. Featuring a patented, first-to-market locking mechanism to encourage patient compliance and anti-microbial Sorbatex padding, Miami J Select is designed to support successful patient outcomes.

## FEATURES AND BENEFITS

- ① Easy height adjustability with a tamper-proof locking mechanism – **Left to Lock**
- ② Anti-microbial & biocompatible Sorbatex™ padding reduces bacterial count by over 99%
- ③ 7 precise fitting sizes with opportunity to fine-tune
- ④ One-piece, anatomically inspired collar back with Flex Edge® designed to reduce occipital pressure
- ⑤ Easy to fit, don and doff with patented Reproducible Fit option\*
- ⑥ Simplified comfort adjustments with integrated and patented Sternal-Relief Dial
- ⑦ Large tracheal opening – improved function for ICU patients

SUGGESTED  
HCPCS CODE  
L0174



**Blue is You™** – intuitive blue Patient Touch Points and multi-colored padding – send patients home with confidence

## ORDERING INFORMATION

Description	Size	Part Numbers
Miami J Select Collar Set (contains the Miami J Select collar and a set of replacement pads)	One-size adjustable (Universal)	MJSR-101



Frontline Medical Services  
Toll Free: 866-681-0444  
www.frontlinemedicalmil.com



Serve proud. Live well.™



PROUD DISTRIBUTOR OF