

REBOUND® POST-OP KNEE

Easy to fit. Easy to use.

Orders may be placed through the DLC Remote Order Entry System (ROES) online ordering application. For assistance, contact DLC Customer Service Support at 303-273-6200.

Contract No: 36C791-20-D-0009



The Össur Logo, trademarks and images are used with permission of Össur hf.



Frontline Medical Services
Toll Free: 866-681-0444
www.frontlinemedicalmil.com



Serve proud. Live well.™

PROUD DISTRIBUTOR OF



REBOUND® POST-OP KNEE

Rebound Post-Op Knee is designed with two priorities in mind: ease of fitting for practitioners and ease of use for patients.

FEATURES AND BENEFITS

- 1 Large buckles featuring hook closures and grip-friendly ridges for simple brace application and removal
- 2 Numbered frame for patient ease-of-use
- 3 Color-coded drop lock and buckles for easy patient education and product use – Blue is You™
- 4 Lock-out ROM pull tabs offer -10° to 120° of ROM control, with OneTouch™ drop-lock from -10° to 30°
- 5 Anti-Migration System wrap for secure suspension
- 6 Quick-fit buckle with a lever arm helps ensure optimal strap tension and brace suspension
- 7 Telescoping frame features easy-to-engage OneTouch™ button with secure locking feature
- 8 Alligator strap tabs enable easy strap trimming
- 9 Moldable, anatomically contoured struts help ensure a precise fit

SMART-FIT™ PACKAGING

Rebound Post-op Knee features Smart-Fit™ packaging to aid in proper brace positioning and fit. Simply open the box, lay it flat, slide it underneath the patient's leg and follow the fitting instructions printed on the box.

ORDERING INFORMATION

Size	Thigh Circumference*	Telescoping	Cool
Universal	16"-32" (41-81cm)	18.5"-27.5" (47-69.5cm)	B-8381480

*Measurement taken at 6" (15cm) above mid-patella



Frontline Medical Services
Toll Free: 866-681-0444
www.frontlinemedicalmil.com



Serve proud. Live well.™



PROUD DISTRIBUTOR OF