

Available through VA Denver Logistics Center (DLC)

Contract No: 36C79120D0009

Orders may be placed through the DLC Remote Order Entry System (ROES) online ordering application. For assistance, contact DLC Customer Service Support at 303-273-6200.

X-Wool Prosthetic Sock with X-STATIC® Pure Silver Fiber

This sock features X-STATIC® – The Silver Fiber®, a high-performance fiber system made with pure silver. Silver is naturally antimicrobial, enabling X-STATIC to inhibit odor in the garment. Silver also has thermodynamic properties, so patients stay cooler in the summer and warmer in winter.

FEATURES

- Thick, absorbent virgin wool provides cushioning
- Stretchy polyester/X-STATIC®/Lycra® core-spun yarn for better fit
- Contours to the limb without stretching out
- Fleeced lining adds softness
- X-STATIC® silver fibers assist in preventing the growth of odor causing bacteria and fungi
- X-STATIC® silver fibers are moisture-wicking and thermodynamic, keeping you cooler in the summer and warmer in the winter
- Available in 3-Ply and 5-Ply thickness
- Machine washable and resists shrinkage

INDICATIONS

- Used to maintain fit of transtibial and transfemoral prostheses

CONTRAINDICATIONS

If rash, blisters, or skin irritations occur, consult your physician and discontinue use of the sock.



KNIT-RITE®
INCORPORATED



Frontline Medical Services
Toll Free: 866-681-0444
www.frontlinemedicalmil.com



Serve proud. Live well.™

X-WOOL® | PROSTHETIC SOCK



X-WOOL® PROSTHETIC SOCK SIZING GUIDE REGULAR

3-Ply Regular	5-Ply Regular	Size	Fits Sock Size	Length
1WX3NAXS	1WX5NAXS	Narrow X-Short	00, 0, AS, AL	8" - 10" (20-25cm)
1WX3NASH	1WX5NASH	Narrow Short	00, 0, AS, AL	10" - 12" (25-30cm)
1WX3NARG	1WX5NARG	Narrow Regular	00, 0, AS, AL	14" - 16" (35-40cm)
1WX3NALG	1WX5NALG	Narrow Long	00, 0, AS, AL	18" - 20" (45-50cm)
1WX3MDXS	1WX5MDXS	Medium X-Short	1, 2	8" - 10" (20-25cm)
1WX3MDSH	1WX5MDSH	Medium Short	1, 2	10" - 12" (25-30cm)
1WX3MDRG	1WX5MDRG	Medium Regular	1, 2	14" - 16" (35-40cm)
1WX3MDLG	1WX5MDLG	Medium Long	1, 2	18" - 20" (45-50cm)
1WX3MDXL	1WX5MDXL	Medium X-Long	1, 2	22" - 24" (55-60cm)
1WX3MDXX	1WX5MDXX	Medium XX-Long	1, 2	26" - 28" (65-70cm)
1WX3WDSH	1WX5WDSH	Wide Short	3, 4	10" - 12" (25-30cm)
1WX3WDMD	1WX5WDMD	Wide Medium	3, 4	14" - 16" (35-40cm)
1WX3WDLG	1WX5WDLG	Wide Long	3, 4	18" - 20" (45-50cm)
1WX3WDXL	1WX5WDXL	Wide X-Long	3, 4	22" - 24" (55-60cm)
1WX3WDXX	1WX5WDXX	Wide XX-Long	3, 4	26" - 28" (65-70cm)
1WX3EWSH	1WX5EWSH	X-Wide Short	5, 6	12" - 14" (30-35cm)
1WX3EWMD	1WX5EWMD	X-Wide Medium	5, 6	16" - 18" (40-45cm)

Thickness Color Coding: 3-Ply = Yellow, 5-Ply = Green
 HCPC: L8420



Place orders directly through the VA Denver Logistics Center (DLC).
 Contact DLC Customer Service Support at 303-273-6200.



Frontline Medical Services
 Toll Free: 866-681-0444
www.frontlinemedicalmil.com



Serve proud. Live well.™