

# X-WOOL® | PROSTHETIC SOCK - HOLE IN TOE

Available through VA Denver Logistics Center (DLC)

Contract No: 36C79120D0009

Orders may be placed through the DLC Remote Order Entry System (ROES) online ordering application. For assistance, contact DLC Customer Service Support at 303-273-6200.

## X-Wool Prosthetic Sock with Hole In Toe for Suspension

*This sock features X-STATIC® – The Silver Fiber®, a high-performance fiber system made with pure silver. Silver is naturally antimicrobial, enabling X-STATIC to inhibit odor in the garment. The hole in toe for suspension is reinforced to prevent fraying.*

### FEATURES

- Reinforced hole in toe to accommodate pin on locking liner
- Thick, absorbent virgin wool provides cushioning
- Stretchy polyester/X-STATIC®/Lycra® core-spun yarn for better fit
- Contours to the limb without stretching out
- Fleeed lining adds softness
- X-STATIC® silver fibers assist in preventing the growth of odor causing bacteria and fungi
- X-STATIC® silver fibers are moisture-wicking and thermodynamic, keeping you cooler in the summer and warmer in the winter
- Available in 3-Ply and 5-Ply thickness
- Machine washable and resists shrinkage

### INDICATIONS

- Maintain fit of transtibial and transfemoral prostheses
- Accommodate for pin on locking liner

### CONTRAINDICATIONS

If rash, blisters, or skin irritations occur, consult your physician and discontinue use of the sock.



**KNIT-RITE®**  
INCORPORATED



Frontline Medical Services  
Toll Free: 866-681-0444  
[www.frontlinemedicalmil.com](http://www.frontlinemedicalmil.com)



Serve proud. Live well.™



## X-WOOL® PROSTHETIC SOCK SIZING GUIDE HOLE-IN-TOE

3-Ply Hole In Toe	5-Ply Hole In Toe	Size	Fits Sock Size	Length
1WX3SNXS	1WX5SNXS	Narrow X-Short	00, 0, AS, AL	8" - 10" (20-25cm)
1WX3SNSH	1WX5SNSH	Narrow Short	00, 0, AS, AL	10" - 12" (25-30cm)
1WX3SNRG	1WX5SNRG	Narrow Regular	00, 0, AS, AL	14" - 16" (35-40cm)
1WX3SNLG	1WX5SNLG	Narrow Long	00, 0, AS, AL	18" - 20" (45-50cm)
1WX3SMXS	1WX5SMXS	Medium X-Short	1, 2	8" - 10" (20-25cm)
1WX3SM SH	1WX5SM SH	Medium Short	1, 2	10" - 12" (25-30cm)
1WX3SMRG	1WX5SMRG	Medium Regular	1, 2	14" - 16" (35-40cm)
1WX3SMLG	1WX5SMLG	Medium Long	1, 2	18" - 20" (45-50cm)
1WX3SWSH	1WX5SWSH	Wide Short	3, 4	10" - 12" (25-30cm)
1WX3SWMD	1WX5SWMD	Wide Medium	3, 4	14" - 16" (35-40cm)

Thickness Color Coding: 3-Ply = Yellow, 5-Ply = Green

HCPC: L8420



**KNIT-RITE®**  
INCORPORATED

Place orders directly through the VA Denver Logistics Center (DLC).  
Contact DLC Customer Service Support at 303-273-6200.



Frontline Medical Services  
Toll Free: 866-681-0444  
[www.frontlinemedicalmil.com](http://www.frontlinemedicalmil.com)

